# **ARSIM2CARE APPLICATION**

## **USER'S MANUAL**

---

**Verónica Rita Dias Coutinho**  
**Paula Escalada-Hernández**  
**Maria da Conceição Giestas Baía Saraiva**  
**Luís Leitão Sarnadas**  
**Rui Carlos Negrão Baptista**  
**Klaas de Proost**  
**Isabel San Martín-Erice**  
**Leticia San Martín-Rodríguez**  
**Nelia Soto-Ruiz**  
**Sandra Tricas-Sauras**  
**Filip Wybraeke**

---

Manual developed as intellectual output from the Project ARsim2care funded by the Program Erasmus + of the European Union.

<http://www.ARsim2care.com/>

The Project is developed in collaboration with three European high education institutions: Escola Superior de Enfermagem de Coimbra (ESEnfC), Erasmushogeschool Brussel (EhB) and Public University of Navarre (UPNA); and the company Industrial Augmented Reality (iAR).

You will find further information about this project and all the materials needed for the correct use of the ARsim2care application, on the Erasmus + website:

<https://ec.europa.eu/programmes/erasmus-plus/projects/eplus-project-details/#project/2017-1-ES01-KA203-038514>

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## ARSIM2CARE APPLICATION USER'S MANUAL

Text © The Editors 2020

First Edition 2020

ISBN 978-84-09-28859-5

This work is published under a Creative Commons Attribution-NonCommercial-NoDerivatives 4.0 International (CC BY-NC-ND 4.0) (<https://creativecommons.org/licenses/by-nc-nd/4.0/>). This license allows to share, copy and redistribute the material for non-commercial purposes and give credit to the original authors and the source, providing a link to the Creative Commons license.

License: CC BY-NC-ND 4.0

---

**Suggested citation:** Dias Coutinho, V. R., Escalada-Hernández, P., Giestas Baía Saraiva, M.G., Leitão Sarnadas, L., Negrão Baptista, R. C., de Proost, K., San Martín-Erice, I., San Martín-Rodríguez, L., Soto-Ruiz, N., Tricas-Sauras, S., & Wybraeke, F. (2020). Arsim2care Application User's Manual. Available at <https://hdl.handle.net/2454/39170>. Licensed under a Creative Commons Attribution-NonCommercial-NoDerivatives 4.0 International (CC BY-NC-ND 4.0).

## CONTRIBUTORS

**Verónica Rita Dias Coutinho RN, PhD**

Adjunct Professor, Medical-Surgical Nursing  
Nursing School of Coimbra  
Coimbra, PORTUGAL

**Paula Escalada Hernández RN, MSc, PhD**

Associate professor  
Department of Health Sciences  
Public University of Navarre  
Pamplona, Navarra, SPAIN

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Pamplona, Navarra, SPAIN

**Sandra Tricas-Sauras, RN, MSc, MRes, PhD**

International Academic Coordinator

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## English

### TABLE OF CONTENTS

ARSIM2CARE APPLICATION: USER'S MANUAL.....	8
APP OPENING ICON.....	8
HOME SCREEN.....	9
MANUAL PROCEDURES SCREEN.....	11
SIMULATION PROCEDURES SCREEN .....	12
RECOGNITION MODE SCREEN .....	13
MANUAL MODE SCREEN .....	17
TEACHING PROTOCOL FOR THE ADEQUATE USE OF THE ARSIM2CARE APP.....	20
PRELIMINARY KNOWLEDGE: INTRODUCTION TO HOOLENS GLASSES USE.....	20
1- Quick capsule on Augmented Reality .....	20
2- About Microsoft HoloLens glasses:.....	22
3- About the ARsim2care Application .....	24
SIMULATION SESSION: PERFORMING A PROCEDURE.....	25
1- Procedure performance with ARsim2care Application.....	25
2- Debriefing.....	26
TECHNICAL ANNEX 1 .....	27
TECHNICAL ANNEX 2 .....	28
TECHNICAL ANNEX 3 .....	29

## ENGLISH

# ARSIM2CARE APPLICATION: USER'S MANUAL

ARsim2care allows the integration of Augmented Reality (AR) with Microsoft HoloLens glasses for clinical procedures simulation, such as: Intramuscular injection, Endotracheal Intubation, Nasogastric Tube Insertion and Suctioning via a Tracheostomy Tube. The Application (ARsim2care App) also provides guidance on the previously mentioned procedures. The use of the application is very intuitive.

### APP OPENING ICON

After Windows is opened in your HoloLens glasses, the ARsim2care App must be opened through the following icon (Fig. 1):

Figure 1. The app opening icon.

**Important:** for a proper functioning of the App, your HoloLens should be linked to a WI-FI connection.

## HOME SCREEN

You will find all the procedures included in the App on the home screen. It appears in the centre of the image, occupying most of the screen (Fig. 2).


Figure 2. The home screen with the procedures included.

At the upper right corner of the screen you will find two icons:

Settings: 1) change language 2) Bluetooth (BT) name: refers to the BT network to which the HoloLens will be linked (in case there are more than one mannequin sensorized) (Fig. 3). 3) BT Debug button: should only be active when the support team requests it (Fig. 4).

Figure 3. The options Setting


Figure 4. Bluetooth Debug.

Access to Procedures Manual.

## MANUAL PROCEDURES SCREEN

Once you access the procedures manual, our App shows you a specific screen in which the navigation will guide you to the different sections and procedures (Fi. 5).


Figure 5. Manual procedures screen.

On the right side of the screen, the icons will guide you through the different sections of the program:

- Definitions
- Indications
- Contraindications
- Materials
- Steps
- Complications

The user can navigate, go forward or backwards through the different sections. This can be done by using the left /right arrows keys or using the voice commands: "NEXT" and "BACK".

### SIMULATION PROCEDURES SCREEN

Each one of the procedure's screens has a similar structure. If the procedure can be performed with a specific anatomic model or with a mannequin, the App gives the user the possibility to choose among either of them (Fig. 6).


Figure 6. Simulation procedures screen.

In the lower part of the screen, there are two icons: the user can choose between them (recognition or manual modes) by moving the central key from one side to another.

“Recognition mode”: The App will work by recognizing the anatomic model and projecting the 3D images over it.

“Manual Mode”: Application will project the 3D images outside the model so the user can manually place and fix them wherever he or she considers.

Central button to switch between one mode or another.

**Important:** for a correct functioning of the glasses' recognition:

- Moderate ambient light: avoid too much or too low light
- Try to create color contrast between the mannequin and the support's surface. It is recommended to use a green or blue surgical cloth under the mannequin to enhance this contrast.
- Recognition angle: it is important to get into a correct position. In any case, this position should be standing up at the feet of the mannequin looking towards its head.

### RECOGNITION MODE SCREEN

While choosing Recognition Mode, the HoloLens will project an anatomic silhouette on your mannequin. This should fit with the anatomical model or the simulation mannequin. You should be careful during this procedure (Fig.7).


Figure 7. Recognition mode screen.

Consequently, HoloLeens will recognize the model and will project the 3D images on the right mannequin section (Fig. 8).



Figure 8. The 3D images on the right mannequin section.

When the images are projected, the user should fix them using the voice command “STOP” so that they fit over the model (or mannequin). After launching the command, the user has to wait for 2-3 seconds before moving so that the image will not move.

Once you get into “Recognition Mode” a floating quick access menu is also displayed (Fig. 9).

Figure 9. The floating quick access menu.

This round button has a central icon (a “house”) for a direct exit to the home page, and four more icons around it (Fig. 10).

Figure 10. Menu button.

Access to the Procedures Manual.

Bluetooth on/off connection (it is recommended to have it disconnected if sensors are not being used).

Manual mode.

Recognition by QR code: the marker must be placed at the anatomic area where the procedure will be developed. It should have the same size as the QR code attached at the Annex 3 of the Manual (6x6cms). The same QR code can be used for all the procedures (Fig. 11) .


Figure 11. QR code.

In case we get into recognition by QR code, the floating quick access menu button changes slightly. This new button gives the opportunity to go again to “Recognition Mode” or to “Manual Mode” (Fig. 12).


Figure 12. Menu button.

### MANUAL MODE SCREEN

While accessing “Manual Mode”, HoloLens will project 3D images on the air, which you will see near you. The user can move those images with hand movements, changing them from one place to another, making them smaller or bigger, as his or her own need (Fig. 13) .

- To turn the image: use the squares located at the central area of the cube edges.
- To make the image bigger or smaller: use the circles located at the cube edges.


Figure 13. The virtual images with the cube edges.

In this case the floating quick access menu button shows the following icons (Fig. 14):


Figure 14. Menu button.

 Access to the Procedures Manual

 Bluetooth on/off connection

 Recognition mode

 Manual mode. Removes interaction cube

**Important:** For the correct functioning of the App it is important to close it once you have finished using it by doing the “bloom gesture”, and clicking the house icon of the menu. In case the screen becomes white, click the x at the upper right part of the screen (Fig. 15).


Figure 15. The finish screen.

# TEACHING PROTOCOL FOR THE ADEQUATE USE OF THE ARSIM2CARE APP

ARsim2care App allows the integration of Augmented Reality (AR) with Microsoft HoloLens glasses for clinical procedures simulation such as: Intramuscular injection, Endotracheal Intubation, Nasogastric Tube Insertion and Suctioning via Tracheostomy Tube.

Using this application implies a radical change in the way of using low fidelity simulation for teaching the above-mentioned procedures. For this reason, it is necessary to establish a teaching protocol which describes the best way to approach a simulation session while using the ARsim2care App to conduct specific procedures. The aim of this protocol is to ensure a most satisfactory learning experience.

## **PRELIMINARY KNOWLEDGE: INTRODUCTION TO HOOLENS GLASSES USE**

Before the beginning of a simulation session with ARsim2care App, it is important that students understand AR technology and know the basic functioning of the HoloLens glasses, as well as the use of our ARsim2care App. For those reasons, a 15 minutes long introduction session is proposed. The group size could be as big as is needed or as big as the classroom capacity allows.

*Preliminary Session for ARsim2care use: ***Introduction to HoloLens glasses use.****

- *Course length: 15 min*
- *Material Resources: Power Point projection*
- *Group size: unlimited*
- *Training content:*

## **1- Quick capsule on Augmented Reality**

In a very simple way, the concept of Augmented Reality can be compared to the one of Virtual Reality (VR). As VR implies the total immersion into a virtual world with a pair of glasses which isolate the user from the real world, AR includes virtual images into our world. This means that using AR, the user does not stop seeing the real world but an amplified one, where certain images appear with the goal to “enhance” the information that the real world already provides.

Some specific examples will help the student understand this technology, for instance, the popular game “Pokemon Go”; but also other applications which are commonly used nowadays based on this technology. For instance, smartphone Apps that super-position over the users’ face virtual elements (i.e. dog’s ears, lips, glasses, etc.), and so on. Those examples will help remembering that AR always works with the help of a device like a smartphone, a tablet, or HoloLens (Fig. 16).


Figure 16. Virtual reality and augmented reality.

Another important concept to share with students is AR’s mechanism of action. Virtual images appear once the device used detects a QR code, or recognizes a specific image or an object: the sooner the student realizes this, the better he will understand the software. The software, as it will be further clarified later, starts working at the very moment that it recognizes the anatomical model over which simulation will be performed. In the case that this recognition is not working, it is also possible to use the App by QR code markers (Fig. 17).


Figure 17. Some examples.

## 2- About Microsoft HoloLens glasses:

Beyond the technology incorporated in this device (that is not necessary to go deeper into it), it is convenient to explain some user's instructions.

First, it is important to explain how to place the HoloLens on the head. Specifically, HoloLens 1, (we are aware that HoloLenses 2 are already in the market, but have a slightly different clamping mechanism). It is also important to know that the glasses' weight is sustained by a headband which allows glasses to fit better. The headband must be put on a proper way: from forehead to nape (as close as possible to the neck) (Fig. 18).

Figure 18. The correct place

It is important to inform students that the **visual field** of the glasses is narrower, and smaller than ours, so if the student remains very close to the virtual image, they will not be able to see it completely.

Additionally, it is important to remember that HoloLens work either by **voice commands** or by **hand gestures**. As per voice commands, they simply need to know the ones relevant for the ARsim2care App (they will be shown further on).

The hand gestures that must be used are worth explaining further in detail, to make the usage of the HoloLens glasses easier and enjoyable. It is thus important to inform the student that:

1. HoloLens detect our hands only when they are in front of our field of vision.
2. There are some specific gestures to perform different actions:
  - To open or close the different applications (“bloom gesture”) (Fig. 19).
  - To activate specific elements use the “mouse pointer” (as clicking with the mouse), as we know it, is a white point that will appear in front of our eyes, and it will shift while moving the head. For clicking it, the hand should be placed in front of the visual field, index finger should be lifted up and lowered down quickly (when glasses detect the finger, the circle increases its size) (Fig. 20).
  - To drag items: using the thumb, index and the middle finger making the gesture of “grabbing something” and “dropping something” (Fig. 21).


Figure 19. The bloom gesture 1

Figure 20. The mouse pointer

Figure 21. The gesture of dropping something.

### 3- About the ARsim2care Application

Finally, it is convenient to show how the basic structure of the App works, so that students are already acquainted with its functioning while using it during the different procedures.

First, how to locate the ARsim2care app once the Windows HoloLens menu is opened. You will find the ARsim2care app in the list provided (alphabetically), as shown in the picture (Fig. 22):


Figure 22. The app opening icon.

## **SIMULATION SESSION: PERFORMING A PROCEDURE**

Simulation sessions must be structured depending on the number of the attended students and depending on the number of available glasses . As an example, the session could be scheduled as follows:

### ***Simulation Session: Performing a procedure.***

- *Course length: 60 minutes.*
- *Resources: HoloLens glasses and simulation material.*
- *Group size: 5-10 students (depending on the procedure).*
- *Training content:*

From a pedagogical point of view, it seems more interesting for students to perform the full procedure with the help of HoloLens, through the ARsim2care App. However, if it is a complex procedure, or if the students' experience with the HoloLens is not optimal, we would suggest to first conduct the procedure without the HoloLens in order to run the session more dynamically.

### **1- Procedure performance with ARsim2care Application**

The sequence of the procedure completion by the student will run as follows:

1. Selection of procedure's material
2. Arrangement of the material and its preparation for the procedure.
3. Placing properly the HoloLens.
4. Placing sterile gloves in case the procedure requires it.
5. Opening the ARsim2care App and selecting the procedure accordingly.
6. Moving forward through the menu until 3D images are settled and fixed at the right place of projection (mannequin/air).
7. Selecting the accurate material to perform the procedure.
8. Performing the procedure.

Once finished, it is compulsory to disinfect the glasses by each student following the protocol anti COVID-19, as established by each educational centre.

Students not performing the procedure will be observing their colleagues, and completing a checklist to evaluate the correct developing of the technique.

## 2- Debriefing

Reflection group:

- Steps performed properly.
- Steps performed incorrectly.
- Difficulties found.
- Advantages found with the use of the glasses.

## TECHNICAL ANNEX 1

### SOFTWARE INSTALLING

- Turn on HoloLens glasses.
- Go to *Settings >> Network & Internet >> WiFi* : Connect HoloLens to the same WiFi network as the computer is connected.
- Once connection is settled, select again the same WiFi network as you're already connected and get into: *Advanced Options*.
- Scroll down until you see the IPv4 address and write it down.
- Open a new window at the computer with Chrome browser and write down the IP address at Chrome URL bar.
- If the browser is warning that it is not a private connection (it only happens during the first connection), you get into advanced configuration and you choose to go to the address anyway.
- You will be asked for HoloLens's User and Password: in case you don't have it or you don't remember it, you can reinstaure it by failing three consecutive times. HoloLens glasses will reset and you can establish a new user.
- To upload the App to the HoloLens glasses, once you're at *Windows Device Portal*, you choose *View >> Apps >> Deploy Apps*, and select: *Choose file*.
- Select: *ARsim2care.appx*
- Wait until installation is finished.
- You already have the ARsim2care application available in your HoloLens glasses.
- To take into account:
  - It is the same procedure for a Software update.
  - It's not convenient being at *Device Portal* and using the ARsim2care app at the same time. It slows the functioning of the App

## TECHNICAL ANNEX 2

### SENSORIZATION

For having a proper visualization of the inner structures while performing the procedures, sensorization of the mannequin's digestive and respiratory tubes are needed.



In the following link you'll find a tutorial video for the correct mannequin sensorization:  
<https://youtu.be/Oa9VI3IVwv0>

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