SUPPLEMENTARY MATERIAL

Effective Degradation of Methylene Blue in Aqueous Solution Using Pd-Supported Cu-doped Ti-Pillared Montmorillonite Catalyst

Annu Joseph^{a,}, Kannan Vellayan^{a,*}, Beatriz González^b, Miguel A. Vicente^{b,*}, Antonio Gil

E-mail addresses: kannanpvl@gmail.com (K. Vellayan), mavicente@usal.es (M.A. Vicente).

^a Department of Chemistry, Government College Kattappana, Kerala, India-685508

^b GIR-QUESCAT, Departamento de Química Inorgánica, Universidad de Salamanca, 37008 Salamanca, Spain

^c INAMAT-Departamento de Ciencias, Universidad Pública de Navarra, 31006 Pamplona, Spain

^{*} Corresponding authors.




Fig S1. UV-Vis absorption spectra of MB solution treated with different doses of MTiCuPd500 catalyst, in the presence of NaBH₄, for 30 minutes.




Fig S2. UV-Vis absorption spectra of MB solution treated with 5 mg of MTiCuPd500 catalyst, in the presence of NaBH₄, for different times.




Fig S3. Evolution with time of UV-Vis absorption spectra of MB solution in the presence of NaBH₄ but without addition of catalyst.




Fig S4. Evolution with time of UV-Vis absorption spectra of MB solution in the presence of MTiCuPd500 catalyst but in the absence of NaBH₄.




Fig. S5. Evolution with time of UV-Vis absorption spectra of MB solution catalyzed by MTiCuPd500 (5 mg) in the presence of NaBH₄ but in absence of light. For illustrative purposes, the spectra in this Figure are plotted in the 200-800 nm range, although the peak in the UV region was not used for calibration of the intensity-concentration relationship.




Fig. S6. Evolution with time of UV-Vis absorption spectra of MB solution catalyzed by MTi500 catalyst (5 mg) in the presence of NaBH₄.